

Animal Surgical Instruments

Intellibio distributes RWD surgical instruments cover surgical scissors, surgical forceps, orthopedic instruments, scalpels and blades, vascular forceps, needle forceps, with more than 1 000 kinds of models.

Surgery Instrument Kit ►

SP0001-G General Surgery Instrument Kit for Rodents

Cat.No.	Product Description	Qty
S11001-08	Spring Scissors (Triangular)-S/S Str/5*0.1mm/8.5cm	1
S11002-08	Spring Scissors (Triangular)-S/S Cvd/5*0.1mm/8.5cm	1
S12009-11	IRIS-Fine Scissors (Round Type)-S/S Str/11.5cm	1
S12010-11	IRIS-Fine Scissors (Round Type)-S/S Cvd/11.5cm	1
S12003-09	IRIS-Fine Scissors (Round Type)-S/S Str/9.5cm	1
S12004-09	IRIS-Fine Scissors (Round Type)-S/S Cvd/9.5cm	1
F11001-11	Micro Forceps-Str, 0.2mm Tips, 11cm	1
F11004-11	Micro Forceps-Cvd, 0.2mm Tips, 11.5cm	1
F12010-10	Dressing Forceps-Str, 1.9mm Tips, 10.5cm	1
F12011-10	Dressing Forceps-Cvd, 1.9mm Tips, 10.5cm	1
F21001-12	HALSTED Artery Forceps-Str, 2.0mm Tips, 12.5cm	2
F21002-12	HALSTED Artery Forceps-Cvd, 2.0mm Tips, 12.5cm	2
F22002-10	HARTMAN Mosquito Forceps-Str, 1.0mm Tips, 10cm	2
F22003-10	HARTMAN Mosquito Forceps-Cvd, 1.0mm Tips, 10cm	2
F31047-12	OLSEN-HEGAR Needle Holders with Scissors-Str, 12cm	1
S13029-14	MAYO-STILLE Dissecting & Operating Scissors-B/B Str/14cm	1
S15001-11	SPENCER Ligature Scissors (Slender Type)-11cm	1
S21010-15	BOHLER Rongeurs (Double)-Light Cvd/2mm Cup/15.5cm	1
S22004-11	Bone Cutters with Flat Blades (SGL)-11.5cm	1
S14014-12	Operating Scissors (Round Type)-S/S Str/12.5cm	1
S31011-01	11# Scalpel Blades (Box of 100pcs)	1
S32001-12	Scalpel Handles 3# Solid-12cm	1
S33009-06	Scalpel Blade Remover-6.5cm	1
F13019-12	ADSON 1x2 Teeth Tissue Forceps, 1.5mm, 12cm	1
R21027-12	STEVENS Hooks, 1 Angled Tooth (5mm long), 12.5cm	1
R22009-01	ALM 4x4 Teeth Retractors-Blunt, 7cm	1
R22005-45	3x3 Teeth Retractors-Blunt, 4.5cm	1
R51003-11	Micro Spatula w/tip W. 2mm x Thk. 0.3mm - Cvd/11cm	1
F35205-60	Sutures w/Needle-#3/8/2.5x7/30cm/6-0	10
F35305-50	PGA Sutures w/Needle-#1/2/4x10/90cm/5-0	10
R31006-06	SS Micro Clamps-Cvd/L*W 6*1mm/17mm	5
R31005-06	SS Micro Clamps-Str/L*W 6*1mm/15mm	5
R34001-14	Clip Applicator for R31005- and R31006-Clamps-14cm	1
SP0000-P	Instrument Storage Portfolio, 31*21cm	1



SP0003-R Microsurgery Instrument Kit for Rat, Cat, Rabbit etc.

Cat.No.	Product Description	Qty
S12009-11	IRIS-Fine Scissors (Round Type)-S/S Str/11.5cm	1
S12010-11	IRIS-Fine Scissors (Round Type)-S/S Cvd/11.5cm	1
F31047-12	OLSEN-HEGAR Needle Holders with Scissors-Str, 12cm	1
F11001-11	Micro Forceps-Str, 0.2mm Tips, 11cm	1
R31005-06	SS Micro Clamps-Str/L*W 6*1mm/15mm	5
R34001-14	Clip Applicator for R31005- and R31006-Clamps-14cm	1
S11001-08	Spring Scissors (Triangular)-S/S Str/5*0.1mm/8.5cm	1
S11002-08	Spring Scissors (Triangular)-S/S Cvd/5*0.1mm/8.5cm	1
F22006-12	HALSTED Mosquito Forceps-Str, 1.0mm Tips, 12.5cm	1
F22007-12	HALSTED Mosquito Forceps-Cvd, 1.0mm Tips, 12.5cm	1
F35401-50	Needle-#3/8/3x10/90cm/5-0 (50/Box)	10
S32001-12	Scalpel Handles 3# Solid-12cm	1
S31011-01	11# Scalpel Blades (Box of 100pcs)	1
SP0000-P	Instrument Storage Portfolio, 31*21cm	1



SP0003-M Microsurgery Instrument Kit for Mice

Cat.No.	Product Description	Qty
S12003-09	IRIS-Fine Scissors (Round Type)-S/S Str/9.5cm	1
S12004-09	IRIS-Fine Scissors (Round Type)-S/S Cvd/9.5cm	1
F31047-12	OLSEN-HEGAR Needle Holders with Scissors-Str, 12cm	1
F11001-11	Micro Forceps-Str, 0.2mm Tips, 11cm	1
R31005-04	SS Micro Clamps-Str/L*W 4*0.75mm/13mm	5
R34001-14	Clip Applicator for R31005- and R31006-Clamps-14cm	1
S11001-08	Spring Scissors (Triangular)-S/S Str/5*0.1mm/8.5cm	1
S11002-08	Spring Scissors (Triangular)-S/S Cvd/5*0.1mm/8.5cm	1
F22002-10	HARTMAN Mosquito Forceps-Str, 1.0mm Tips, 10cm	1
F22003-10	HARTMAN Mosquito Forceps-Cvd, 1.0mm Tips, 10cm	1
F35401-50	Needle-#3/8/3x10/90cm/5-0 (50/Box)	10
S32001-12	Scalpel Handles 3# Solid-12cm	1
S31011-01	11# Scalpel Blades (Box of 100pcs)	1
SP0000-P	Instrument Storage Portfolio, 31*21cm	1



SP0004-M
Brain Microdialysis Surgery Instrument Kit for Mice

Cat.No.	Product Description	Qty
S31011-01	11# Scalpel Blades (Box of 100pcs)	1
S32001-12	Scalpel Handles 3# Solid-12cm	1
S33007-12	ZIEGLER Knives-4mm Cutting	1
S12003-09	IRIS-Fine Scissors (Round Type)-S/S Str/9.5cm	1
S12004-09	IRIS-Fine Scissors (Round Type)-S/S Cvd/9.5cm	1
S14014-11	Operating Scissors (Round Type)-S/S Str/11.5cm	1
F12005-10	IRIS Dissecting Forceps-Str, 0.8mm Tips, 10cm	2
F12006-10	Forceps-Light Cvd, 0.8mm Tips, 10cm	1
F22002-10	HARTMAN Mosquito Forceps-Str, 1.0mm Tips, 10cm	2
F31047-12	Needle Holders with Scissors-Str, 12cm	1
F35205-60	Sutures w/Needle- ϕ 3/8/2.5 $\times$ 7/30 mm /6-0 (50/Box)	10
F22003-10	HARTMAN Mosquito Forceps-Cvd, 1.0mm Tips, 10cm	2
SP0000-P	Instrument Storage Portfolio, 31 $\times$ 21cm	1


SP0004-R
Brain Microdialysis Surgery Instrument Kit for Rat

Cat.No.	Product Description	Qty
S31011-01	11# Scalpel Blades (Box of 100pcs)	1
S32001-12	Scalpel Handles 3# Solid-12cm	1
S33007-12	ZIEGLER Knives-4mm Cutting	1
S12005-10	IRIS-Fine Scissors (Round Type)-S/S Str/10.5cm	1
S12006-10	IRIS-Fine Scissors (Round Type)-S/S Cvd/10.5cm	1
S14014-12	Operating Scissors (Round Type)-S/S Str/12.5cm	1
F12010-10	Dressing Forceps-Str, 1.9mm Tips, 10.5cm	2
F12011-10	Dressing Forceps-Cvd, 1.9mm Tips, 10.5cm	1
F22006-12	HALSTED Mosquito Forceps-Str, 1.0mm Tips, 12.5cm	2
F31047-12	OLSEN-HEGAR Needle Holders with Scissors-Str, 12cm	1
F35205-60	w/Needle- ϕ 3/8/2.5 $\times$ 7/30 mm /6-0 (50/Box)	10
R22009-01	ALM 4x4 Teeth Retractors-Blunt, 7cm	1
F22007-12	HALSTED Mosquito Forceps-Cvd, 1.0mm Tips, 12.5cm	2
SP0000-P	Instrument Storage Portfolio, 31 $\times$ 21cm	1


SP0002-G
General Ophthalmic Surgery Instrument Kit for Rodents

Cat.No.	Product Description	Qty
S11001-08	Spring Scissors (Triangular)-S/S Str/5 $\times$ 0.1mm/8.5cm	1
S11002-08	Spring Scissors (Triangular)-S/S Cvd/5 $\times$ 0.1mm/8.5cm	1
S11003-08	Spring Scissors (Triangular)-S/S Angled/5 $\times$ 0.1mm/8.5cm	1
S12005-10	IRIS-Fine Scissors (Round Type)-S/S Str/10.5cm	1
S12006-10	IRIS-Fine Scissors (Round Type)-S/S Cvd/10.5cm	1
S15006-09	O'BRIEN Stitch Scissors-S/S Angled/9.5cm	1
F12013-10	Dressing Forceps w/o Serrations-45 $^{\circ}$ Cvd, S/S, 10cm	1
F13019-12	ADSON 1x2 Teeth Tissue Forceps, 1.5mm, 12cm	1
F12005-10	IRIS Dissecting Forceps-Str, 0.8mm Tips, 10cm	1
F12006-10	IRIS Dissecting Forceps-Light Cvd, 0.8mm Tips, 10cm	1
F18011-09	Forceps w/5mm Platform-Str, 0.4mm Tips, 9cm	1
F31047-12	HEGAR Needle Holders with Scissors-Str, 12cm	1
F35401-50	Needle- ϕ 3/8/3 $\times$ 10/90 mm /5-0 (50/Box)	10
R21027-12	STEVENS Hooks, 1 Angled Tooth (5mm long), 12.5cm	1
R22029-04	Colibri Eye Specula-1.5cm spread, 4cm	1
SP0000-P	Instrument Storage Portfolio, 31 $\times$ 21cm	1


SP0007-M
Universal dissecting Surgery Instrument Kit for Mice

Cat.No.	Product Description	Qty
S12003-09	IRIS-Fine Scissors (Round Type)-S/S Str/9.5cm	1
S12004-09	IRIS-Fine Scissors (Round Type)-S/S Cvd/9.5cm	1
S13003-11	Fine Dissecting Scissors (Flat Type)-B/B Str/11.5cm	1
S13004-11	Fine Dissecting Scissors (Flat Type)-B/B Cvd/11.5cm	1
S15001-09	SPENCER Ligature Scissors (Slender Type)-9cm	1
S32001-12	Scalpel Handles 3# Solid-12cm	1
S31015-01	15# Scalpel Blades (Box of 100pcs)	1
F12005-10	IRIS Dissecting Forceps-Str, 0.8mm Tips, 10cm	2
F12006-10	IRIS Dissecting Forceps-Light Cvd, 0.8mm Tips, 10cm	2
F22002-10	HARTMAN Mosquito Forceps-Str, 1.0mm Tips, 10cm	2
F22003-10	HARTMAN Mosquito Forceps-Cvd, 1.0mm Tips, 10cm	2
F13029-10	IRIS 1x2 Teeth Tissue Forceps-Str, 0.8mm Tips, 10cm	1
F31047-12	OLSEN-HEGAR Needle Holders with Scissors-Str, 12cm	1
F35305-50	PGA Sutures w/Needle- ϕ 1/2/4 $\times$ 10/90 mm /5-0 (50/Box)	10
R22005-45	3x3 Teeth Retractors-Blunt, 4.5cm	1
R31005-04	SS Micro Clamps-Str/L*W 4 $\times$ 0.75mm/13mm	5
R42002-12	Spinal Cord Hook - Tip Dia. 3mm/12cm	1
R34001-14	Clip Applicator for R31005- and R31006-Clamps-14cm	1
SP0000-P	Instrument Storage Portfolio, 31 $\times$ 21cm	1


SP0005-R
Brain Microinjection Surgery Instrument Kit for Rat

Cat.No.	Product Description	Qty
S31011-01	11# Scalpel Blades (Box of 100pcs)	1
S32001-12	Scalpel Handles 3# Solid-12cm	1
S12005-10	IRIS-Fine Scissors (Round Type)-S/S Str/10.5cm	1
S12006-10	IRIS-Fine Scissors (Round Type)-S/S Cvd/10.5cm	1
F12010-10	Dressing Forceps-Str, 1.9mm Tips, 10.5cm	2
F12011-10	Dressing Forceps-Cvd, 1.9mm Tips, 10.5cm	2
F22006-12	HALSTED Mosquito Forceps-Str, 1.0mm Tips, 12.5cm	2
F22007-12	HALSTED Mosquito Forceps-Cvd, 1.0mm Tips, 12.5cm	2
F31047-12	OLSEN-HEGAR Needle Holders with Scissors-Str, 12cm	1
S14014-12	Operating Scissors (Round Type)-S/S Str/12.5cm	1
F35205-60	Sutures w/Needle- ϕ 3/8/2.5 $\times$ 7/30 mm /6-0 (50/Box)	10
R22009-01	ALM 4x4 Teeth Retractors-Blunt, 7cm	1
R51003-11	Micro Spatula w/tip W. 2mm x Thk. 0.3mm - Cvd/11cm	1
SP0000-P	Instrument Storage Portfolio, 31 $\times$ 21cm	1


SP0005-M
Brain Microinjection Surgery Instrument Kit for Mice

Cat.No.	Product Description	Qty
S31011-01	11# Scalpel Blades (Box of 100pcs)	1
S32001-12	Scalpel Handles 3# Solid-12cm	1
S12003-09	IRIS-Fine Scissors (Round Type)-S/S Str/9.5cm	1
S12004-09	IRIS-Fine Scissors (Round Type)-S/S Cvd/9.5cm	1
S14014-11	Operating Scissors (Round Type)-S/S Str/11.5cm	1
F12005-10	IRIS Dissecting Forceps-Str, 0.8mm Tips, 10cm	2
F12006-10	IRIS Dissecting Forceps-Light Cvd, 0.8mm Tips, 10cm	2
F22002-10	HARTMAN Mosquito Forceps-Str, 1.0mm Tips, 10cm	2
F31047-12	OLSEN-HEGAR Needle Holders with Scissors-Str, 12cm	1
F35205-60	Sutures w/Needle- ϕ 3/8/2.5 $\times$ 7/30 mm /6-0 (50/Box)	10
F22003-10	HARTMAN Mosquito Forceps-Cvd, 1.0mm Tips, 10cm	2
R51003-11	Spatula w/tip W. 2mm x Thk. 0.3mm - Cvd/11cm	1
SP0000-P	Instrument Storage Portfolio, 31 $\times$ 21cm	1



SP0006-M
Bone Surgery Instrument Kit for Mice

Cat.No.	Product Description	Qty
S14014-11	Operating Scissors (Round Type)-S/S Str/11.5cm	1
S16014-09	NAIL Scissors (Broad Type)-S/S Str/9cm	1
S21020-14	Friedman-Pearson Rongeurs (SGL)-Str/0.7mm Cup/14cm	1
S22004-11	Bone Cutters with Flat Blades (SGL)-11.5cm	1
S21023-14	Friedman-Pearson Rongeurs (SGL)-Cvd/0.7mm Cup/14cm	1
S23007-12	LAMBOTTE Osteotomes - 4mm Cutting Edge/12.5cm	1
S33006-13	GRAEFE Scalpels-22mm Cutting Edge/13cm	1
SP0000-P	Instrument Storage Portfolio, 31*21cm	1



SP0006-R
Bone Surgery Instrument Kit for Rat, Cat, Rabbit etc.

Cat.No.	Product Description	Qty
S14014-14	Operating Scissors (Round Type)-S/S Str/14cm	1
S16014-09	NAIL Scissors (Broad Type)-S/S Str/9cm	1
S21021-14	Pearson Rongeurs (SGL)-Str/1.0mm Cup/14cm	1
S22005-12	Bone Cutters with Flat Blades (SGL)-12.5cm	1
S23007-12	LAMBOTTE Osteotomes - 4mm Cutting Edge/12.5cm	1
S23002-16	FREER Chisels - Cvd/4mm Cutting Edge/16cm	1
S33006-13	GRAEFE Scalpels-22mm Cutting Edge/13cm	1
S21024-14	Pearson Rongeurs (SGL)-Cvd/1.0mm Cup/14cm	1
SP0000-P	Instrument Storage Portfolio, 31*21cm	1



SP0007-C
Univeral Dissecting Surgery Instrument Kit for Cat, Rabbit

Cat.No.	Product Description	Qty
S12009-11	RIS-Fine Scissors (Round Type)-S/S Str/11.5cm	1
S12010-11	IRIS-Fine Scissors (Round Type)-S/S Cvd/11.5cm	1
S12031-11	Standard Scissors (Round Type)-S/S Full Cvd/11.5cm	1
S13005-14	Dissecting Scissors (Round Type)-B/B Str/14.5cm	1
S13006-14	Dissecting Scissors (Round Type)-B/B Cvd/14.5cm	1
F12010-13	Dressing Forceps-Str, 2.1mm Tips, 13cm	2
F12011-13	Dressing Forceps-Cvd, 2.1mm Tips, 13cm	2
F21001-12	HALSTED Artery Forceps-Str, 2.0mm Tips, 12.5cm	2
F21002-12	HALSTED Artery Forceps-Cvd, 2.0mm Tips, 12.5cm	2
S15001-11	SPENCER Ligature Scissors (Slender Type)-11cm	1
S32001-12	Scalpel Handles 3# Solid-12cm	1
S31015-01	15# Scalpel Blades (Box of 100pcs)	1
F13014-12	1x2 Teeth Tissue Forceps-Str, 1.0mm Tips, 12.5cm	1
F31047-14	OLSEN-HEGAR Needle Holders with Scissors-Str, 14cm	1
F35305-50	PGA Sutures w/Needle-1/2/4x10/90cm/5-0 (50/Box)	10
R22015-02	WULLSTEIN 3x3 Teeth Retractors-Blunt, 13cm	1
R31008-26	Schwartz Micro Clamps-Str/26mm	5
R42002-12	Spinal Cord Hook - Tip Dia. 3mm/12cm	1
SP0000-P	Instrument Storage Portfolio, 31*21cm	1



SP0007-R
Univeral Dissecting Surgery Instrument Kit for Rat

Cat.No.	Product Description	Qty
S12005-10	IRIS-Fine Scissors (Round Type)-S/S Str/10.5cm	1
S12006-10	IRIS-Fine Scissors (Round Type)-S/S Cvd/10.5cm	1
S12010-11	IRIS-Fine Scissors (Round Type)-S/S Cvd/11.5cm	1
S13005-14	Dissecting Scissors (Round Type)-B/B Str/14.5cm	1
S13006-14	Dissecting Scissors (Round Type)-B/B Cvd/14.5cm	1
F21001-12	HALSTED Artery Forceps-Str, 2.0mm Tips, 12.5cm	2
F21002-12	HALSTED Artery Forceps-Cvd, 2.0mm Tips, 12.5cm	2
F12010-10	Dressing Forceps-Str, 1.9mm Tips, 10.5cm	2
F12011-10	Dressing Forceps-Cvd, 1.9mm Tips, 10.5cm	2
F13029-10	IRIS 1x2 Teeth Tissue Forceps-Str, 0.8mm Tips, 10cm	1
F31047-12	OLSEN-HEGAR Needle Holders with Scissors-Str, 12cm	1
S15001-11	SPENCER Ligature Scissors (Slender Type)-11cm	1
S32001-12	Scalpel Handles 3# Solid-12cm	1
S31015-01	15# Scalpel Blades (Box of 100pcs)	1
F35305-50	PGA Sutures w/Needle-1/2/4x10/90cm/5-0 (50/Box)	10
R22009-01	ALM 4x4 Teeth Retractors-Blunt, 7cm	1
R31005-06	SS Micro Clamps-Str/L*W 6*1mm/15mm	5
R42002-12	Spinal Cord Hook - Tip Dia. 3mm/12cm	1
R34001-14	Applicator for R31005- and R31006-Clamps-14cm	1
SP0000-P	Instrument Storage Portfolio, 31*21cm	1



SP0008-M
Cardio-cerebral Perfusion Surgery Instrument Kit for Mice

Cat.No.	Product Description	Qty
S14014-11	Operating Scissors (Round Type)-S/S Str/11.5cm	1
S12003-09	IRIS-Fine Scissors (Round Type)-S/S Str/9.5cm	1
S11001-08	Spring Scissors (Triangular)-S/S Str/5*0.1mm/8.5cm	1
F12007-10	IRIS Dissecting Forceps-Large Cvd, 0.8mm Tips, 10cm	1
F12010-09	Dressing Forceps-Str, 1.8mm Tips, 9cm	2
F22002-10	HARTMAN Mosquito Forceps-Str, 1.0mm Tips, 10cm	2
F22003-10	HARTMAN Mosquito Forceps-Cvd, 1.0mm Tips, 10cm	2
F35305-50	PGA Sutures w/Needle-1/2/4x10/90cm/5-0 (50/Box)	10
F35205-60	Sutures w/Needle-3/8/2.5x7/30cm/6-0 (50/Box)	10
R21027-12	STEVENS Hooks, 1 Angled Tooth (5mm long), 12.5cm	1
R22005-45	3x3 Teeth Retractors-Blunt, 4.5cm	1
F31047-12	OLSEN-HEGAR Needle Holders with Scissors-Str, 12cm	1
R31005-04	SS Micro Clamps-Str/L*W 4*0.75mm/13mm	5
R34001-14	Clip Applicator for R31005- and R31006-Clamps-14cm	1
R42002-12	Spinal Cord Hook - Tip Dia. 3mm/12cm	1
SP0000-P	Instrument Storage Portfolio, 31*21cm	1



SP0008-R
Cardio-cerebral Perfusion Surgery Instrument Kit for Rat

Cat.No.	Product Description	Qty
S14014-12	Operating Scissors (Round Type)-S/S Str/12.5cm	1
S12005-10	IRIS-Fine Scissors (Round Type)-S/S Str/10.5cm	1
S11002-08	Spring Scissors (Triangular)-S/S Str/5*0.1mm/8.5cm	1
F12007-10	IRIS Dissecting Forceps-Large Cvd, 0.8mm Tips, 10cm	1
F12005-10	IRIS Dissecting Forceps-Str, 0.8mm Tips, 10cm	2
F22006-12	HALSTED Mosquito Forceps-Str, 1.0mm Tips, 12.5cm	2
F22007-12	HALSTED Mosquito Forceps-Cvd, 1.0mm Tips, 12.5cm	2
F35305-50	PGA Sutures w/Needle-1/2/4x10/90cm/5-0 (50/Box)	10
F35205-60	Sutures w/Needle-3/8/2.5x7/30cm/6-0 (50/Box)	10
R21027-12	STEVENS Hooks, 1 Angled Tooth (5mm long), 12.5cm	1
R22005-45	3x3 Teeth Retractors-Blunt, 4.5cm	1
F31047-12	OLSEN-HEGAR Needle Holders with Scissors-Str, 12cm	1
R31005-06	SS Micro Clamps-Str/L*W 4*0.75mm/13mm	5
R42002-12	Spinal Cord Hook - Tip Dia. 3mm/12cm	1
R34001-14	Clip Applicator for R31005- and R31006-Clamps-14cm	1
SP0000-P	Instrument Storage Portfolio, 31*21cm	1



SP0009-R General Surgery Instrument Kit for Rat

Cat.No.	Product Description	Qty
F11020-11	Micro Forceps-Str, 0.05mmx0.01mm Tips, 11cm	1
F11021-11	Micro Forceps-45°Angled, 0.10x0.06mm Tips, 11cm	1
S32001-12	Scalpel Handles 3# Solid-12cm	1
S31011-01	11# Scalpel Blades (Box of 100pcs)	1
F35401-50	Needle-3/8/3×10/90 /5-0 (50/Box)	10
F35205-60	Needle-3/8/2.5×7/30 /6-0 (50/Box)	10
F31047-12	OLSEN-HEGAR Needle Holders with Scissors-Str, 12cm	1
S11001-08	Spring Scissors (Triangular)-S/S Str/5*0.1mm/8.5cm	1
S11002-08	Spring Scissors (Triangular)-S/S Cvd/5*0.1mm/8.5cm	1
R22009-01	ALM 4x4 Teeth Retractors Blunt, 7cm	1
F12005-10	IRIS Dissecting Forceps-Str, 0.8mm Tips, 10cm	1
F12006-10	IRIS Dissecting Forceps-Light Cvd, 0.8mm Tips, 10cm	1
F12007-10	IRIS Dissecting Forceps-Large Cvd, 0.8mm Tips, 10cm	1
F13029-10	IRIS 1x2 Teeth Tissue Forceps-Str, 0.8mm Tips, 10cm	1
S12009-11	IRIS-Fine Scissors (Round Type)-S/S Str/11.5cm	1
S12010-11	IRIS-Fine Scissors (Round Type)-S/S Cvd/11.5cm	1
F21001-12	HALSTED Artery Forceps-Str, 2.0mm Tips, 12.5cm	1
F21002-12	HALSTED Artery Forceps-Cvd, 2.0mm Tips, 12.5cm	1
SP0000-P	Instrument Storage Portfolio, 31*21cm	1



SP0009-M General Surgery Instrument Kit for Mouse

Cat.No.	Product Description	Qty
F11020-11	Micro Forceps-Str, 0.05mmx0.01mm Tips, 11cm	1
F11021-11	Micro Forceps-45°Angled, 0.10x0.06mm Tips, 11cm	1
S32001-12	Scalpel Handles 3# Solid-12cm	1
S31011-01	11# Scalpel Blades (Box of 100pcs)	1
F35401-50	NA Terylene Sutures w/Needle-3/8/3×10/90/5-0 (50/Box)	10
F31047-12	OLSEN-HEGAR Needle Holders with Scissors-Str, 12cm	1
S11001-08	VANNAS Spring Scissors (Triangular)-S/S Str/5*0.1mm/8.5cm	1
S11002-08	VANNAS Spring Scissors (Triangular)-S/S Cvd/5*0.1mm/8.5cm	1
R22029-04	Colibri Eye Specula-1.5cm spread, 4cm	1
F12005-10	IRIS Dissecting Forceps-Str, 0.8mm Tips, 10cm	1
F12006-10	IRIS Dissecting Forceps-Light Cvd, 0.8mm Tips, 10cm	1
F12007-10	IRIS Dissecting Forceps-Large Cvd, 0.8mm Tips, 10cm	1
F13029-10	IRIS 1x2 Teeth Tissue Forceps-Str, 0.8mm Tips, 10cm	1
S12003-09	IRIS-Fine Scissors (Round Type)-S/S Str/9.5cm	1
S12004-09	IRIS-Fine Scissors (Round Type)-S/S Cvd/9.5cm	1
F22002-10	HARTMAN Mosquito Forceps-Str, 1.0mm Tips, 10cm	1
F22003-10	HARTMAN Mosquito Forceps-Cvd, 1.0mm Tips, 10cm	1
SP0000-P	Instrument Storage Portfolio, 31*21cm	1



SP0010-B Basic Surgical Instrument Kit for Rat & Mice

Cat.No.	Product Description	Qty
S11001-08	VANNAS Spring Scissors (Triangular)-S/S Str/5*0.1mm/8.5cm	1
S12005-10	IRIS-Fine Scissors (Round Type)-S/S Str/10.5cm	1
F12005-10	IRIS Dissecting Forceps-Str, 0.8mm Tips, 10cm	1
F12006-10	IRIS Dissecting Forceps-Light Cvd, 0.8mm Tips, 10cm	1
F11020-11	Micro Forceps-Str, 0.05mmx0.01mm Tips, 11cm	1
F22002-10	HARTMAN Mosquito Forceps-Str, 1.0mm Tips, 10cm	1
F22003-10	HARTMAN Mosquito Forceps-Cvd, 1.0mm Tips, 10cm	1
F31047-12	OLSEN-HEGAR Needle Holders with Scissors-Str, 12cm	1
S33007-12	ZIEGLER Knives-4mm Cutting Edge/1.5mm Width/12cm	1
R31005-06	SS Micro Clamps-Str/L*W 6*1mm/15mm	1
R31006-06	SS Micro Clamps-Cvd/L*W 6*1mm/17mm	1
SP0000-P	Instrument Storage Portfolio, 31*21cm	1



sterilizer

In animal experiments, surgical instruments are easy to be contaminated, rapid sterilization of surgical instruments allows the operation to continue, but the traditional sterilization methods (such as high-pressure steam sterilization) sterilization process and cooling time longer, delay and affect the experimental process. Rapid sterilizer can solve this problem for customers. It can kill germs, pathogens and other microorganisms on the surface of surgical instruments in a very short period of time (about 15-20 seconds), so as to achieve the purpose of disinfection and sterilization.

1. Glass bead heating from room temperature to target temperature (240-270 DEG C) time no more than 20 minutes
2. With power indicator light, power 110V and 220v;



Surgical Instruments

Micro Spring Scissors



Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S11001-08	Straight Sharp&Sharp	7.5mm	1.7mm	8.5cm
S11002-08	Curved Sharp&Sharp	6mm	1.4mm	8.5cm
S11003-08	Curved Sharp&Sharp	5mm	1.3mm	8cm



Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S11035-08	Straight Sharp&Sharp	3.5mm	1.3mm	8cm
S11036-08	Straight Sharp&Sharp	9mm	2.0mm	8.5cm
S11037-08	Curved Sharp&Sharp	8mm	1.9mm	8.5cm

Micro Spring Scissors



Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S11005-09	Curved Blunt&Blunt	10mm	2.6mm	9cm



Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S11029-12	Curved Blunt&Blunt	11mm	1.8mm	11cm



S11024-11

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S11024-11	Curved Sharp&Sharp	18mm	3.7mm	11.5cm



S11014-12

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S11014-12	Straight Sharp&Sharp	10mm	3.2mm	12.5cm



S11007-12


Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S11007-12	Straight Sharp&Sharp	20mm	3.8mm	11.5cm

Delicate Scissors



S12003-09

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S12003-09	Straight Sharp&Sharp	25mm	5.0mm	10cm
S12005-10	Straight Sharp&Sharp	22mm	4.8mm	10.5cm
S12009-11	Straight Sharp&Sharp	25mm	5.8mm	11.5cm
S12014-13	Straight Sharp&Sharp	31mm	8mm	13cm



S12004-09

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S12004-09	Curved Sharp&Sharp	20mm	4.9mm	9.5cm
S12006-10	Curved Sharp&Sharp	18mm	4.3mm	10.5cm
S12010-11	Curved Sharp&Sharp	24mm	5mm	11cm
S12015-13	Curved Sharp&Sharp	26mm	7.2mm	13cm




S12011-12



S12012-12

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S12011-12	Straight Sharp&Sharp	31mm	5.6mm	12cm
S12012-12	Curved Sharp&Sharp	31mm	5.6mm	12cm



S12029-11




S12030-11



S12031-11

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S12029-11	Straight Sharp&Sharp	23mm	5mm	11.5cm
S12030-11	Curved Sharp&Sharp	23mm	5.3mm	11.5cm
S12031-11	Curved Sharp&Sharp	28mm	5.9mm	11.5cm
S12034-12	Straight Sharp&Sharp	16mm	5.5mm	12.5cm



S12047-10

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S12047-10	Straight Sharp&Sharp	18mm	4.5mm	10cm



S18001-09



S18002-09

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S18001-09	Straight Sharp&Sharp (Teeth)			9cm
S18001-10	Straight Sharp&Sharp (Teeth)			10.5cm
S18001-11	Straight Sharp&Sharp (Teeth)	25mm	5mm	11.5cm
S18002-09	Curved Sharp&Sharp (Teeth)	20mm	4.9mm	9cm
S18002-10	Curved Sharp&Sharp (Teeth)	22mm	5mm	10.5cm
S18002-11	Curved Sharp&Sharp (Teeth)	23mm	5mm	11.5cm




S18004-09



S18005-09

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S18004-09	Curved Sharp&Sharp	21mm	4.8mm	9cm
S18004-10	Curved Sharp&Sharp	22mm	5mm	10.5cm
S18004-11	Curved Sharp&Sharp	22mm	5.2mm	11.5cm
S18005-09	Straight Sharp&Sharp	20mm	4.5mm	9cm
S18005-10	Straight Sharp&Sharp	22mm	5mm	10.5cm
S18005-11	Straight Sharp&Sharp	22mm	5mm	11.5cm



S13001-10

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S13001-10	Straight Sharp&Sharp	26mm	5.8mm	10cm




S13003-11



S13004-11


Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S13003-11	Straight Blunt&Blunt	25mm	6.4mm	11.5cm
S13004-11	Curved Blunt&Blunt	28mm	6.6mm	11.5cm



S13020-11

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S13020-11	Straight Sharp&Sharp	20mm	5.2mm	12cm

Dissecting Scissors




S13019-12

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S13019-12	Curved Blunt&Blunt	25mm	4.6mm	12.5cm



Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S13027-18	Straight Blunt&Blunt	32mm	7mm	18cm



S13048-14

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S13048-14	Curved Blunt&Blunt	32mm	9mm	14cm



S13035-16

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S13035-16	Straight Blunt&Blunt	45mm	9.5mm	16.5cm



S13047-14

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S13047-14	Straight Blunt&Blunt	40mm	9.8mm	14.5cm



S13029-14

S13030-14

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S13029-14	Straight Blunt&Blunt	39mm	9.2mm	14cm
S13030-14	Curved Blunt&Blunt	41mm	9.4mm	14cm



S13031-17

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S13031-17	Straight Blunt&Blunt	44mm	14mm	17cm



S13038-14

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S13038-14	Straight Blunt&Blunt	41mm	10mm	14.5cm

Surgical Scissors



S14011-13

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S14011-13	Straight Sharp&Sharp	28mm	10mm	13cm




S14016-11

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S14016-11	Curved Sharp&Sharp	34mm	7.5mm	11.5cm



S14002-10

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S14002-10	Curved Sharp&Blunt	30mm	7.2mm	10cm
S14002-13	Curved Sharp&Blunt	42mm	8.7mm	13cm
S14002-14	Curved Sharp&Blunt	47mm	9.8mm	14.5cm



S14001-10

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S14001-10	Straight Sharp & Blunt	30mm	7.3mm	10.5cm
S14001-11	Straight Sharp & Blunt	35mm	8.7mm	11.5cm
S14001-13	Straight Sharp & Blunt			13cm
S14001-14	Straight Sharp & Blunt	48mm	9.8mm	14.5cm
S14001-15	Straight Sharp & Blunt	48mm	9.8mm	15cm
S14001-16	Straight Sharp & Blunt	53mm	10mm	16.5cm



S14014-10

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S14014-10	Straight Sharp&Sharp	34mm	7.2mm	10.5cm
S14014-11	Straight Sharp&Sharp	35mm	8.2mm	11.5cm
S14014-12	Straight Sharp&Sharp	38mm	8mm	12.5cm
S14014-13	Straight Sharp&Sharp	42mm	8mm	13cm
S14014-14	Straight Sharp&Sharp	43mm	8.9mm	14cm
S14014-15	Straight Sharp&Sharp	48mm	10mm	15.5cm
S14014-16	Straight Sharp&Sharp	52mm	10mm	16.5cm
S14014-17	Straight Sharp&Sharp	55mm	10.7mm	17.5cm



S14015-10

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S14015-10	Straight Blunt&Blunt	33mm	7.5mm	12cm
S14015-11	Straight Blunt&Blunt	36mm	8.2mm	11.5cm
S14015-12	Straight Blunt&Blunt	37mm	8.7mm	12.5cm
S14015-14	Straight Blunt&Blunt	46mm	9.75mm	14.5cm
S14015-16	Straight Blunt&Blunt	54mm	11.3mm	16.5cm
S14015-17	Straight Blunt&Blunt	53mm	12mm	17.5cm



S14017-10

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S14017-10	Curved Blunt&Blunt	33mm	7.5mm	10.5cm
S14017-16	Curved Blunt&Blunt	51mm	10.3mm	16cm



S14003-14

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S14003-14	Straight sharp & blunt	47mm	9.7mm	14.5cm



S15006-09

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S15006-09	Curved sharp & sharp	23mm		9.5cm

Ligature Scissors



Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S15001-09		22mm	5mm	9cm
S15001-11				11cm

Bone Cutters



S16001-10

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S16001-10	Straight	2mm	9.9mm	10.5cm



S22005-12

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S22005-12		11.5mm	20mm	13cm
S22004-11	plane	16mm		11cm

Bone Chisel



S23002-16



S23007-12

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S23002-16			4mm	16cm
S23007-12			4mm	12.5cm

scalpel




S33007-12



S33006-13

Cat No.	Description of Blade	Cutting Edge (Length)	Cutting Edge (Width)	Total Length
S33007-12		1.5mm	6mm	12cm
S33006-13		25mm	25mm	13cm

Surgical Handles



S32001-12 (matches 10#, 11#, 12#, 15# blades)



S32003-12 (matches 10#, 11#, 12#, 15# blades)



S32004-13 (matches 10#, 11#, 12#, 15# blades)



S32005-12

(matches 10#, 11#, 12#, 15# blades)

Cat No.		Total Length
S32001-12	3#(Standard)	12cm
S32003-12	3#(Ruler)	12cm

Cat No.		Total Length
S32004-13	4#(Standard)	13cm
S32005-12	7#	12.5cm

Surgical Blades



Cat No.		
S31010-01	10#	Box of 100pcs



Cat No.		
S31011-01	11#	Box of 100pcs



Cat No.		
S31012-01	12#	Box of 100pcs



Cat No.		Total Length
S31015-01	15#	Box of 100pcs



Cat No.		
S31023-02	23#	Box of 100pcs



Cat No.		Total Length
S33009-06	Blade Remover	6.5cm

Delicate Forceps



F11001-11

Cat No.	Tip Shape	Tip Length	Tip Width	Tip Thickness	Total Length
F11001-11	Straight		0.20mm		11cm
F11002-11	Straight		0.15mm		11cm
F11001-13	Straight		0.15mm		13cm
F11020-11	Straight	28mm	0.05mm	0.1mm	11cm



F11003-11

Cat No.	Tip Shape	Tip Length	Tip Width	Tip Thickness	Total Length
F11003-11	45°Curved	8mm	0.3mm		11cm
F11021-11	45°Curved	8mm	0.10mm	0.06mm	11cm
F11003-13	45°Curved	8mm	0.15mm		13cm



F13002-11

Cat No.	Tip Shape	Tip Length	Tip Width	Tip Thickness	Total Length
F13002-11	Straight		0.30mm		11cm



F13004-11

Cat No.	Tip Shape	Tip Width	Tip Thickness	Total Length
F13004-11	Curved	0.3mm		11.5cm
F11004-11	Curved	0.2mm		11.5cm
F11022-11	Curved	0.07mm	0.03mm	11cm




F11008-09

Cat No.	Tip Shape	Tip Width	Tip Thickness	Total Length
F11008-09	1x2Teeth 5mm Platform	0.28mm		9.5cm



F11014-11




F11019-11




F11016-11

Cat No.	Tip Shape	Tip Width	Tip Thickness	Total Length
F11014-11	Straight	0.15mm		11cm
F11019-11	Straight	0.3mm		11cm
F11016-11	Straight	0.3mm		11cm



F11027-13



F11028-13

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F11027-13	Straight	26mm		13cm
F11028-13	Curved	32mm		13cm



F11029-11



F11030-11



F11031-11

Cat No.	Tip Shape	Tip Width	Tip Thickness	Total Length
F11029-11	Straight Titanium Alloy	0.05mm	0.01mm	11cm
F11030-11	45°Curved Titanium Alloy	0.10mm	0.06mm	11cm
F11031-11	Curved Titanium Alloy	0.07mm	0.03mm	11cm

Dissecting Forceps



F12002-12



F12003-12

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F12002-12	Dissecting Forceps	13mm	0.6mm	12cm
F12003-12	Dissecting Forceps	18mm	1.5mm	12cm



F13019-12

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F13019-12	Straight 1x2Teeth		1.8mm	12cm



F13027-12

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F13027-12	7x7Teeth Straight		3.0mm	12cm



F12012-10



F12013-10

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F12012-10	Straight Not Teeth		0.4mm	10cm
F12013-10	45°Curved Not Teeth	7mm	0.5mm	10cm



F15007-15

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F15007-15	50°Curved	15mm	1.2mm	15cm



F14009-10



F15017-14

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F14009-10	Curved		8mmx3.3mm	10.5cm
F15017-14	Straight	18mm	2.0mm	14.5cm



F12010-09

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F12010-09	Straight	12mm	1.8mm	9.5cm
F12010-10	Straight	15mm	1.8mm	10.5cm
F12010-11	Straight	16mm	2.5mm	11.5cm
F12010-13	Straight	17mm	1.3mm	13cm



F12011-09

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F12011-09	Curved	11mm	1.8mm	9.5cm
F12011-10	Curved	18mm	1.9mm	10.5cm
F12011-13	Curved	16mm	1.2mm	13cm



F12019-10

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F12019-11	Straight	16.5mm	1.9mm	11.5cm
F12019-13	Straight	17.5mm	2.0mm	13cm
F12019-16	Straight	19mm	12.1mm	16cm



F13029-10



F13030-10

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F13029-10	Straight		1.3mm	10.5cm
F13030-10	Curved		1.2mm	10.5cm
F13009-11	Straight		1.1mm	9cm



F12018-10

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F12018-14	Straight	19mm	3.3mm	14.5cm



F12016-15



F12029-09

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F12016-15	Straight	19mm	1.4mm	15.5cm
F12029-09	Straight	12mm	0.8mm	9.5cm



F13014-12

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F13014-12	Straight 1x2Teeth		1.0mm	13cm



F12005-10



F12006-10



F12007-10



F14012-10

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F12005-10	Straight	13mm	1.0mm	10cm
F12006-10	Curved	12mm	1.0mm	10cm
F12007-10	Curved	10mm	0.9mm	10cm
F14012-10	Curved		0.5mm	10cm

Tissue Forceps



F13023-10

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F13023-14	Straight 1x2Teeth		1.8mm	14cm



F13023-15

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F13028-15	Straight circinal	9.5mm	7.0mm	15cm



F18003-12

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F18003-12	Straight		0.5mm	12cm



F18011-09

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F18011-09	Straight	16mm	0.4mm	9cm

Glass Cover Forceps




F18015-10



F18001-10

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F18015-10	Straight	11mm	4.0mm	10.5cm
F18001-10	Curved	12mm	4.7mm	10.5cm

Tissue Forceps




F18006-13

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F18008-15	Straight 4x5Teeth	39mm	5.0mm	15.0cm
F18009-15	Straight 5x6Teeth	44mm	5.9mm	15.5mm




F21002-12



F21001-12


Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F21001-12	Straight	21mm	1.0mm	12.5cm

Haemostatic Forceps




F21011-14

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F21011-14	Straight	36mm	2.0mm	14.5cm



F21016-14

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F21016-14	Curved	35mm	2.0mm	14.5cm



F21020-16


F21019-16

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F21019-16	Straight	37mm	1.5mm	16cm
F21020-16	Curved	37mm	1.5mm	16cm




F22006-12

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F22006-12	Straight	20mm	1.2mm	12.5cm



F21007-12

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F22007-12	Curved	21mm	1.5mm	12.5cm



F22003-09

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F22003-09	Curved	19mm	1.0mm	9cm
F22003-10	Curved	18mm	1.0mm	10cm



F21004-12




F21003-12

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F21003-12	Straight 1x2Teeth	22mm	1.5mm	12.5cm
F21004-12	Curved 1x2Teeth	20mm	1.6mm	12.5cm




F21006-14




F21005-14

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F21005-14	Straight	35mm	1.5mm	14.5cm
F21006-14	Curved	22mm	1.7mm	14.5cm




F22002-09

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F22002-09	Straight	19mm	1.0mm	9cm
F22002-10	Straight	19mm	1.4mm	10.5cm



F21015-16

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F21015-16	Straight	37mm	1.8mm	16cm



F23022-16

F23021-16

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F23021-16	Straight	25mm	1.7mm	16cm
F23022-16	Curved	36mm	1.5mm	16cm



F31047-12

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F31047-12	Straight	8mm	2.0mm	12cm
F31047-14	Straight	13mm	2.1mm	14cm

Needle Holders



F31007-13




F31002-13

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F31002-13	Straight	14mm	1.0mm	13.5cm
F31007-13	Straight	13mm	1.2mm	13.0cm



F31031-15

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F31031-15	Straight	19mm	2.6mm	15.5cm



F31025-13

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F31025-13	Straight	16mm	2.0mm	13cm




F31034-12

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
F31034-12	Straight	16mm	2.0mm	13cm
F31034-14	Straight	16mm	2.0mm	14cm

suture



Sutures w/Needle- Δ 3/8/2.5 $\times$ 7/30cm/
6-0 (50/Box)
No. F35205-60



PGA Sutures w/Needle- $\circ$ 1/2/4 $\times$ 10/90cm/
5-0 (50/Box)
No. F35305-50




Needle- $\circ$ 3/8/3 $\times$ 10/90cm/5-0 (50/Box)
No. F35401-50

Bone Rongeurs



S21002-15

Cat No.	Jaw Shape	Jaw Size (Cup)	Jaw Length	Total Length
S21002-15	Curved	3.5mm	11mm	15cm



S21010-15

Cat No.	Jaw Shape	Jaw Size (Cup)	Jaw Length	Total Length
S21010-15	Curved	2mm	18mm	15.5cm



S21027-17

Cat No.	Jaw Shape	Jaw Size (Cup)	Jaw Length	Total Length
S21027-17	Straight	1.0mm	12mm	17cm



S21019-14

Cat No.	Jaw Shape	Jaw Size (Cup)	Jaw Length	Total Length
S21019-14	Straight	0.5mm	7mm	14cm
S21020-14	Straight	0.7mm	6.1mm	14cm
S21021-14	Straight	1.0mm	9.5mm	14cm




S21022-14

Cat No.	Jaw Shape	Jaw Size (Cup)	Jaw Length	Total Length
S21023-14	Curved	0.7mm	7mm	14cm
S21024-14	Curved	1.0mm	9mm	14cm

Blood Vessel Clamps




Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
R31005-04	Straight	4mm	0.75mm	16mm
R31005-06	Straight	6mm	1.0mm	19mm
R31005-08	Straight	8mm	2.0mm	17mm
R31005-10	Straight	10mm	2.0mm	20mm
R31006-04	Curved	6.3mm	1.0mm	16mm
R31006-06	Curved	6.6mm	1.0mm	17mm
R33001-38	Straight	11mm	3.5mm	38mm
R33001-48	Straight	14mm	3.7mm	48mm
R33002-38	Curved	11mm	3.75mm	38mm
R33002-48	Curved	15mm	3.7mm	48mm
R33003-01	Straight	17mm	3.7mm	58mm
R33003-02	Curved	18.5mm	4.0mm	58mm
R31008-26	Straight	11mm	2.0mm	28mm
R31009-26	30°Curved	14mm	2.4mm	28mm
R31010-26	45°Curved	15mm	2.5mm	29mm




R34001-14
Used only to clamp R31005, R31006 series vascular clamps

Cat No.	Tip Shape	Tip Length	Tip Width	Total Length
R34001-14		8.4mm		14mm

Retraeters



R22010-01



R22009-01

Cat No.	Teeth Shape	Teeth Length	Teeth Width	Teeth Depth	Maximum Spread	Total Length
R22009-01	Blunt 4x4Teeth	9.5mm	1mm	4.5mm	60mm	7.5cm
R22010-01	Blunt 4x4Teeth	11.5mm	1mm	5.5mm	80mm	10cm




Cat No.	Teeth Shape	Teeth Length	Teeth Width	Teeth Depth	Maximum Spread	Total Length
R21027-12		5mm	1.5mm	5mm		12.5cm

R22029-03



R22012-45




Cat No.	Teeth Shape	Teeth Length	Teeth Width	Teeth Depth	Maximum Spread	Total Length
R22029-03		12mm	0.9mm	6mm	22mm	3cm
R22029-04		12.5mm		6.6mm	29mm	4cm
R22012-45		13mm		7.5mm	32mm	4cm


R22001-50



R22005-45



R22006-55




Cat No.	Teeth Shape	Teeth Length	Teeth Width	Teeth Depth	Maximum Spread	Total Length
R22001-50	Blunt 1x1Teeth	7.5mm	1.5mm	8.9mm	40mm	5cm
R22005-45	Blunt 3x3Teeth	11mm	1mm	7mm	36mm	5cm
R22006-55	Blunt 3x3Teeth	12mm	1mm	8mm	35.5mm	5.5cm




R22011-50

Cat No.	Teeth Shape	Teeth Length	Teeth Width	Teeth Depth	Maximum Spread	Total Length
R22011-50	Sharp 3x3Teeth	10.8mm		8mm	32mm	5.5cm
R22030-03	Sharp 3x3Teeth	4.5mm		2.2mm	21mm	3.5cm



R22013-10

Cat No.	Teeth Shape	Teeth Length	Teeth Width	Teeth Depth	Maximum Spread	Total Length
R22013-11	Blunt 2x3Teeth		13.5x7.8mm	13mm	61mm	11.5cm



R22015-02

Cat No.	Teeth Shape	Teeth Length	Teeth Width	Teeth Depth	Maximum Spread	Total Length
R22015-02	Blunt 3x3Teeth	13.5mm	13mm	0.5mm	62mm	13cm