



## Red Inactinic Desk Lamp



**Desk lamp with red LEDs providing red inactinic lighting, ideal for working in reverse light cycles. The flexible arm allows you to direct the light exactly where it is needed.**

Applications: Work in reverse light cycles (circadian rhythms), observation or laboratory handling.

- Helps preserve the **biological rhythm** of subjects in dark environments
- **Reduces eye strain** and impact on **sleep** when returning to daylight
- Used in chronobiology to study the effects of light on the sleep-wake cycle

### Features

- Red color (narrow spectrum centered at 630 nm)
- Adjustable flexible arm
- Height: 30 cm
- Supplied with a power adapter

### Model

#### Reference

A-1901-00109

#### Description

Red inactinic desk lamp